
# Axio Observer materials Inverted microscope

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Knowledge of this manual is required for the operation of the instrument. Would you therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning safe operation of the instrument. The specifications are subject to change; the manual is not covered by an update service. Unless expressly authorized, forwarding and duplication of this document, as well as utilization and communication of its contents are not permitted. Violations will entail an obligation to pay compensation. All rights reserved in the event of granting of patents or registration of a utility model.

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Carl Zeiss Microscopy GmbH  
Königsallee 9-21  
37081 Göttingen, Germany

Number of this manual: 431004-8044-010  
Date of issue: Version 2, 09/28/2015

**CONTENTS**

	Page
<b>1</b>	<b>GENERAL INFORMATION..... 2</b>
<b>2</b>	<b>SAFETY INSTRUCTIONS ..... 3</b>
2.1	General information ..... 3
2.2	Warranty information ..... 7
2.3	Warning and information labels on the device..... 8
<b>3</b>	<b>DESIGNATED USE ..... 9</b>
<b>4</b>	<b>INSTALLATION NOTES AND FIRST-TIME SET-UP..... 10</b>
4.1	Installation notes ..... 10
4.1.1	Space requirements ..... 10
4.2	First-time set-up ..... 10
4.3	Operating and function controls ..... 11
4.3.1	Axio Observer 3 materials (manual with coded objective nosepiece) ..... 11
4.3.2	Axio Observer 5 materials (encoded, semi-motorized) ..... 13
4.3.3	Axio Observer 7 materials (motorized) ..... 15
4.4	Switching on the instrument ..... 17
4.5	Connectors on the rear of the Axio Observer 3 materials ..... 17
4.6	Connectors on the rear of the Axio Observer 5 materials ..... 18
4.7	Connectors on the rear of the Axio Observer 7 materials ..... 18
4.8	Switching the microscope and the ballast unit for the HBO 100 on and off ..... 19
4.9	Switching off the instrument ..... 19
4.10	Replacing fuses ..... 20
4.10.1	Changing the fuses in the Axio Observer 3 materials and 5 materials ..... 20
4.10.2	Changing the fuses in the Axio Observer 7 materials ..... 21
4.10.3	Changing the fuses in the ballast unit for the HBO 100 ..... 21
<b>5</b>	<b>SHIPPING CONDITIONS..... 22</b>
<b>6</b>	<b>AMBIENT CONDITIONS..... 23</b>
<b>7</b>	<b>TECHNICAL DATA ..... 24</b>

## 1 GENERAL INFORMATION

Axio Observer materials microscopes have been designed, produced and tested in compliance with DIN EN 61010-1 (IEC 61010-1) "Safety requirements for electrical equipment for measurement, control, and laboratory use", meet the requirements of the RoHS Directive 2011/65/EC and are marked with the **CE** mark.

Axio Observer models 3 materials, 5 materials and 7 materials for material applications meet the requirements of the European directives 2006/95/EC, 2014/35/EC "Low Voltage Directives" and 2004/108/EC, 2014/30/EC "EMC directives".

These Safety Instructions contain all information and warnings to be observed by the operator.

The devices shall be disposed of in accordance with the WEEE Directive 2012/19/EC.

The following warning and information symbols are used in this Operating Manual:

**CAUTION**

*This symbol indicates a potential hazard to the user.*

**CAUTION**

*Disconnect the instrument from the power supply before opening!*

**CAUTION**

*Optical radiation is emitted.*

*Do not look into the laser beam! It may be injurious to the eyes.*

**CAUTION**

*Energy-rich UV radiation!*

*Risk of injury to the eyes and skin!*

**CAUTION**

*Inflammable substances, fire hazard!*

**ATTENTION**

*This symbol indicates a potential hazard to the instrument or system.*

**ATTENTION**

*Crushing hazard*

**ATTENTION**

*Read the Operating Manual*

**NOTE**

*This symbol indicates an instruction which requires particular attention.*

**ATTENTION**

*Stand-by: The device is not disconnected from the power supply.*

**ATTENTION**

*Hot surface!*

## 2 SAFETY INSTRUCTIONS

### 2.1 General information

Regulations for occupational health and safety must be observed when operating the HXP 120 V source of UV radiation. The national legal requirements must also be observed.

- Follow the operating manual for the HXP 120 V illuminator as supplied by the manufacturer.
- The lamp must be replaced according to the manufacturer's instructions. There is otherwise a risk of burning or explosion when the lamp is replaced.

- The user should not expose him/herself to UV radiation from the HXP 120 V illuminator.

- Never look directly into the optical fiber when the HXP 120 V illuminator is switched on. Failure to observe this precaution may result in irreversible eye injuries!

Avoid touching the hot lamp housing. Always pull out the power plug before changing the lamps and allow the instrument to cool down for approx. 15 minutes.

The Axio Observer may not be put into operation if the dust cover is still on. Before putting on the dust cover, always check that the device is switched off and has cooled down, otherwise this could cause a fire.

Immersol 518 F<sup>®</sup> immersion oil is a skin irritant. Avoid contact with skin, eyes and clothing.

In the event of contact with the skin, rinse with plenty of water and soap. In the event of contact with the eyes, rinse the eye immediately with plenty of water for at least five minutes. If irritation persists, consult a physician.

Proper disposal of Immersol 518 F<sup>®</sup> immersion oil: Ensure that immersion oil does not enter surface water or the sewage system.

It is essential to read the safety data sheet provided with Immersol 518 F<sup>®</sup> immersion oil.

Gas discharge lamps such as the HBO 100 emit ultraviolet radiation which may cause burns to eyes and skin. Never look directly into the light and avoid direct skin exposure. When using the microscope, always use the device's protective equipment (e.g. special attenuation filters).

Gas discharge lamps develop high internal pressure when hot. They should therefore only be replaced when cool. Protective gloves and eyewear should be worn (for detailed information, please see operating manual 423010-7144-001).

If there are reflector modules in the reflector turret which are fitted with neutral beam splitters or partial mirrors in the beam-splitting mirror position or if these are fitted in the reflector turret itself, looking into the eyepiece can cause injury to the eye if an HBO, X-Cite or HXP lamp is switched on. This applies particularly if the specimen or the specimen holder has reflective properties. Appropriate radiation attenuation measures must be taken (e.g. the use of neutral filters) to prevent damage to the eyes.

When using fluorescence filters, the heat protection filter, which provides protection from radiant heat emitted by the microscope illuminator, must not be removed. Fluorescence filters are heat-sensitive and their function may be impaired if the heat protection filter is removed.

The reflected and transmitted illumination lamps are assigned to LED risk category 2 under DIN EN 62471:2009. Do not look directly into the illumination light.

The desktop power supply units which are available as a microscopy accessory permit mains voltages in the 100 to 240 V  $\pm 10\%$ , 50/60 Hz range without the voltage setting on the instrument having to be changed.

Desktop power supply units available as microscopy accessories should not come into contact with moisture. If the housing is damaged, the desktop power supply unit should be taken out of use. The microscope must only be operated with the desktop power supply unit contained in the scope of delivery.

The operating manuals of the Axio Observer 431004-7144-001, the light sources and the software must be carefully observed.

If it is determined that protective measures are no longer effective, the instrument must be switched off and secured against inadvertent operation. Please contact ZEISS Service or the ZEISS Microscopy Service to repair the instrument.

Please note that the Axio Observer is a precision opto-electronic device. Incorrect handling can easily impair the functioning of the device, or even damage it, and will render any warranty claims invalid.

Switching the Axio Observer 5 materials, 7 materials off using the standby button only turns off the internal computer. The mains supply is not switched off.

The Axio Observer 7 materials can only be disconnected from the mains supply using the on/off switch on the external VP232-2 power unit. Disconnect Axio Observer 5 materials from the power supply.

Keep hands out of the path of the motorized X/Y scanning stages when these are being moved or positioned.

Do not dispose of defective instruments in regular domestic waste; these should be disposed of in accordance with the prevailing legal requirements.

Specimens should also be disposed of in compliance with the prevailing legal requirements and internal operating procedures.

Axio Observer and its original accessories are to be used for the microscopy procedures described in these safety instructions only. The manufacturer cannot assume any liability for any other applications of the instrument, including the use of individual modules or components.

Modifications and repairs to this instrument and any devices operated in combination with the Axio Observer are to be carried out by our service department or by authorized personnel only. The manufacturer accepts no liability for damage caused by unauthorized access to the interior of the instrument. Failure to comply with this shall also render all warranty claims invalid.

The Axio Observer microscopes 3 materials and 5 materials are equipped with a power supply unit integrated in the stand allowing line voltages to be used in the ranges 100 to 127 V and 200 to 240 V AC  $\pm 10\%$ , 50/60 Hz, without the voltage setting on the instrument having to be changed.

Voltage is supplied to the Axio Observer 7 materials by the external VP232-2 power supply unit that belongs to the stand. This is designed for use in the ranges 100 to 240 V  $\pm 10\%$ , 50/60 Hz, without the voltage setting on the instrument having to be changed.

The ballast units for the HBO 100 (ebqb 100-04-z) and XBO 75 (ebx 75-04-z) and HXP 120 are designed for use in the range from 100 to 240 V AC, 50 – 60 Hz. The devices automatically adapt to the mains voltage.

Do not install the device close to heat sources such as radiators, or in direct sunlight. Broad temperature fluctuations and vibrations should be avoided.

The power plugs may only be connected to sockets with earth contact. The protective capacity must not be rendered ineffective by the use of extension cables with no grounding conductor.

Do not replace detachable power cables with power cables that do not meet the specifications. Only the specified power cables should be used.

Before connecting the power cable, ensure that the mains voltage matches that on the type plate of the Axio Observer.

Only specially authorized experts or service staff are permitted to open the device. The Axio Observer may only be operated in enclosed areas.

Disconnect the device from the mains power supply before gaining access to the interior of the device or changing the fuse.

Only use fuses for the given rated current which comply with the specifications on the fuse holders and in this manual. Use of makeshift fuses and short-circuiting of the fuse holders are not permitted.

Do not use portable multiple sockets. Cables are to be laid such that there is no risk of tripping, or are to be covered.

Never disconnect the power cable while the device is in use. Use the power switch to turn the devices off.

To ensure proper function of the device, the Axio Observer should be subjected to an annual safety inspection. The safety checks must be carried out by an authorized ZEISS service technician. All national safety checks should be carried out.

Clogged or covered ventilation slats may lead to heat build-up that will damage the instrument and, in extreme cases, cause a fire. Always keep the ventilation slats clear and ensure that no objects enter or fall into the instrument through the ventilation slats.

Dust and dirt may impair the devices' performance. The devices must be effectively protected from such influences and covered with the dust cover when not in use. This should only be put on once the device has been switched off and has cooled down. Otherwise this could cause a fire.

Do not operate the supplied equipment in areas in which there is a danger of explosion or in the presence of volatile anaesthetics or flammable solvents, such as alcohol, benzene, etc.

Only operate the device on a hard, non-flammable surface.

The devices are not equipped with any special equipment to protect them against corrosive, potentially infectious, toxic, radioactive, or other substances that could be hazardous to health. All legal regulations must be complied with when handling such substances, particularly the prevailing national rules for accident prevention.

To prevent the device from damage, ensure that the generation of electrostatic charges is avoided through suitable design of the workplace environment.

## 2.2 Warranty information

The manufacturer guarantees that the device is free of material or manufacturing defects when delivered. Any defects must be reported immediately and steps taken to minimize damage. If such a defect is reported, the device manufacturer shall be obliged to correct the fault, either by repairing the instrument or replacing it with a new one, at the manufacturer's discretion. No warranty is given for defects caused by natural wear and tear (particularly of wearing parts) and improper use of the device.

The instrument manufacturer shall not be liable for damage caused by misuse, negligence or any other tampering with the instrument, particularly the removal or replacement of instrument components, or the use of accessories from other manufacturers. Such actions shall invalidate any warranty claims.

With the exception of the work described in these safety instructions, no maintenance or repair work is to be carried out on these microscopes. Repairs may only be performed by ZEISS Service or individuals specially authorized by ZEISS Service. In the event of a problem with the instrument, please contact the ZEISS microscopy service team in Germany or your local ZEISS overseas representative.

2.3 Warning and information labels on the device

English


Fig. 1 Warning labels on the microLED attachment lamp



Fig. 2 Labels on the Axio Observer

### 3 DESIGNATED USE

Manufacturer's designation: Inverted reflected light microscope with optional transmitted light equipment

Designation:                   Axio Observer 3 materials (manual / coded version)  
                                      Axio Observer 5 materials (coded / semi-motorized version)  
                                      Axio Observer 7 materials (fully motorized version, incl. motorized Z drive)

Axio Observer materials microscopes 3 materials, 5 materials, 7m are inverted light microscopes for universal use. They are used in all areas of research-based and industrial microscopy. Inverted microscopes permit the unrestricted use of conventional samples as a result of the highly accommodating sample compartment. This facilitates the examination of samples, workpieces etc. of large dimensions.

**Fig. 3** Axio Observer materials microscopes 3 materials, 5 materials and 7 materials

## 4 INSTALLATION NOTES AND FIRST-TIME SET-UP

### 4.1 Installation notes

#### 4.1.1 Space requirements

The following requirements apply regarding installation of the device:

- distance of the system to the wall, at least 3 cm

It must also be ensured that the floor on which the support/table stands is subject to the least possible vibration. Only operate the device on a hard, non-flammable surface.

The parameters described under Ambient conditions and Technical data must also be observed (section 6 and 7).

### 4.2 First-time set-up

Because of the complexity of the equipment and to ensure correct function, the Axio Observer will be installed and set up for use on site by your ZEISS representative.

This includes the following:

- Installation of the microscope, assembly and adjustment of all components (where these are not factory-adjusted)
- Connection of cables and power supply
- Training

Retrofitting and conversion possibilities for the microscope are described in the Axio Observer operating manual 431004-7144-001.

## 4.3 Operating and function controls

### 4.3.1 Axio Observer 3 materials (manual with coded objective nosepiece)


Fig. 4 Operating and function controls of the Axio Observer 3 materials (manual with coded objective nosepiece)

---

**Key to Fig. 4:**

- 1** On / Off switch
- 2** Left sideport
- 3** Coarse / fine focus drive (left side) with finger wheel for fine focus, flat
- 4** Light path switching control (left sideport / vis)
- 5** Microscope stage (with inserted K mounting frame for reflected light and stage pinhole aperture)
- 6** Slot for iris diaphragm slider as reflected light luminous-field diaphragm
- 7** Slot for polarizer slider - A 6x30 mm
- 8** Slot for aperture diaphragm slider MAT or FL attenuator
- 9** Slot for 3-position filter slider, d=25 mm
- 10** Light manager set button
- 11** Drive knobs for controlling XY positioning of the mechanical stage
- 12** 6-position objective nosepiece HD DIC M27 coded
- 13** Compensator mount 6x20
- 14** Reflector turret
- 15** Focus drive coarse / fine (right side)
- 16** RL button for switching the LED or the reflected light shutter (fluorescence) on and off
- 17** Knurled wheel for LED light / halogen illuminator intensity control
- 18** Binocular tube
- 19** Binocular section of the binocular tube
- 20** Eyepiece
- 21** Eyepiece adjustment ring

4.3.2 Axio Observer 5 materials (encoded, semi-motorized)


Fig. 5 Axio Observer 5 materials (encoded, semi-motorized)

**Key to Fig. 5:**

- 1** On / off button
- 2** Left sideport
- 3** Coarse / fine focus drive (left side) with finger wheel for fine focus, flat
- 4** Vertical stop for focus drive
- 5** Light path switching control (left / right sideport / vis)
- 6** Light path switching control (baseport / vis / frontport)
- 7** Microscope stage (with inserted K mounting frame for reflected light and stage pinhole aperture)
- 8** Holder with LCD display
- 9** Slot for iris diaphragm slider as reflected light luminous-field diaphragm (manual or motorized)
- 10** Slot for polarizer slider - A 6x30 mm
- 11** Slot for aperture diaphragm slider MAT or FL attenuator (manual or motorized)
- 12** Slot for 3-position filter slider, d=25 mm
- 13** LM set key
- 14** Drive knobs for controlling XY positioning of the mechanical stage
- 15** Right sideport
- 16** 6-position objective nosepiece HD DIC M27 (coded)
- 17** Compensator mount 6x20
- 18** Reflector turret (coded or motorized)
- 19** Optovar turret control wheel (max. 3 positions)
- 20** Focus drive coarse / fine (right side)
- 21** Control ring, right
- 22** RL button for switching the LED / halogen light or the reflected light shutter (fluorescence) on and off
- 23** Knurled wheel for LED light / halogen illuminator intensity control
- 24** Frontport
- 25** Binocular tube
- 26** Binocular section of the binocular tube
- 27** Eyepiece
- 28** Eyepiece adjustment ring

4.3.3 Axio Observer 7 materials (motorized)


Fig. 6 Axio Observer 7 materials (motorized)

Key to **Fig. 6**:

- 1 Standby button
- 2 Left sideport
- 3 Focus drive coarse / fine (left side)
- 4 Control ring, left
- 5 Microscope stage (with inserted K mounting frame for reflected light and stage pinhole aperture)
- 6 Slot for iris diaphragm slider as reflected light luminous-field diaphragm (manual or motorized)
- 7 Slot for polarizer slider - A 6x30 mm
- 8 Slot for aperture diaphragm slider MAT or FL attenuator (manual or motorized)
- 9 Slot for 3-position filter slider, d=25 mm
- 10 LM set key
- 11 Drive knobs for controlling XY positioning of the mechanical stage
- 12 Right sideport
- 13 6-position objective nosepiece HD DIC M27 mot. ACR (motorized)
- 14 Compensator mount 6x20
- 15 Reflector turret (coded or motorized)
- 16 Coarse / fine focus drive (right side) with finger wheel for fine focus, flat
- 17 Control ring, right
- 18 TFT display
- 19 RL button for switching the LED / halogen light or the reflected light shutter (fluorescence) on and off
- 20 Knurled wheel for LED light / halogen illuminator intensity control
- 21 Frontport
- 22 Binocular tube
- 23 Binocular section of the binocular tube
- 24 Eyepiece
- 25 Eyepiece adjustment ring

A full description of operation is provided in the operating manual of the Axio Observer 431004-7144-001.

#### 4.4 Switching on the instrument

Plug the power cable (supplied) into the power socket of the Axio Observer 3 materials or 5 materials. Switch on the device using the power switch on the microscope,

The Axio Observer 7 materials is supplied with voltage from the external power unit VP232-2.

- Plug connecting plug of power supply VP232-2 into the corresponding connector (Fig. 10/1) on the rear side of the stand.
- Plug the power cable first into the power connector of the power supply VP232-2 and then into a power outlet.
- Turn on the power supply VP232-2 using the power switch (Fig. 7/1).
- Press the standby button on the Axio Observer 7 materials.


Fig. 7 External power unit VP232-2 of the Axio Observer 7 materials

#### 4.5 Connectors on the rear of the Axio Observer 3 materials

Switch off the microscope before connecting any components.

##### Key to figure 8:

- 1 Connector for LED light / halogen light transmitted light (output 1)
- 2 Reflected light/transmitted light selector switch (LED light / halogen light)
- 3 Connector for LED light / halogen light reflected light (output 2)
- 4 Power connector
- 5 CAN connecting socket
- 6 USB port


Fig. 8 Axio Observer 3 materials (rear)


Fig. 9 Axio Observer 5 materials (rear)

#### 4.6 Connectors on the rear of the Axio Observer 5 materials

##### Key to figure 9:

- 1 Connector for LED light / halogen light transmitted light (output 1)
- 2 Reflected light / transmitted light selector switch (LED light / halogen light)
- 3 Connector for LED light / halogen light reflected light (output 2)
- 4 Power socket
- 5 Connector for LCD display  
Transmitted light illumination support
- 6 Connector for transmitted light shutter
- 7 RS-232 port
- 8 TCP/IP port
- 9 USB port
- 10 External high speed shutter connector
- 11 CAN connecting socket
- 12 Trigger socket (IN/OUT) for shutter


Fig. 10 Axio Observer 7 materials (rear)

#### 4.7 Connectors on the rear of the Axio Observer 7 materials

##### Key to figure 10:

- 1 Power supply unit socket VP232-2
- 2 Docking station port (optional)
- 3 Connector for LED light / halogen light transmitted light
- 4 CAN connecting sockets
- 5 not used
- 6 RS-232 port
- 7 Connector for transmitted light shutter
- 8 TCP/IP port
- 9 USB port
- 10 External high speed shutter connector
- 11 CAN connecting socket
- 12 Trigger socket (IN/OUT) for shutter

#### 4.8 Switching the microscope and the ballast unit for the HBO 100 on and off

Axio Observer 3 materials:

- Switch the microscope on and off using the power switch (on the left of the stand, Fig. 11/1).

Axio Observer 5 materials:

- Switch the microscope on and off using the power button (Fig. 11/1).

Axio Observer 7 materials:

- Switch on the microscope's external power supply unit.
- Start up the microscope using the standby button (on the left of the stand, Fig. 11/1).
- To switch the microscope off, press the standby button then switch off the external power supply unit.

The monitoring LED (3, Fig. 11) lights up when the microscope is on.

Ballast unit:

- If a fluorescence illuminator (e.g. HBO 100) is connected, switch the ballast unit on (Fig. 11/2) and off using the power switch.

Further operation descriptions are provided in the operating manual of the Axio Observer 431004-7144-001.


Fig. 11 Axio Observer and ballast unit for HBO 100

#### 4.9 Switching off the instrument

The Axio Observer 3 materials can be disconnected from the mains using the power switch on the microscope.


Use the standby button on the left of the stand of the Axio Observer 5 materials, 7 materials to shut down the internal computer in the microscope. The device is not disconnected from the power supply.

If the Axio Observer 7 materials is not used for a longer period, the external power supply unit 7 materials can be disconnected from the mains using its on/off switch. Disconnect Axio Observer 5 materials from the power supply.

## 4.10 Replacing fuses

Pull out the power plug before changing fuses.

**Fig. 12 Replacing the fuses in the Axio Observer 3 materials and 5 materials**

### 4.10.1 Changing the fuses in the Axio Observer 3 materials and 5 materials

The power fuses for the input voltage are located in the back of the stand of the Axio Observer 3 materials and 5 materials. The fuse compartment is integrated into the power socket and contains two T 5.0 A/H 250 V 5 x 20 mm fuses. If the device fuses fail the cause must first of all be ascertained and technical problems properly remedied.

- Pull out the power plug.
- Pull out the fuse holder (Fig. 12/2) frontward.
- Remove the fuses from the fuse carrier and insert new fuses.
- Reinsert the fuse carrier into the fuse compartment (Fig. 12/1) until it clicks into place.
- Plug in the power cable.

#### 4.10.2 Changing the fuses in the Axio Observer 7 materials

The Axio Observer 7 materials is supplied with voltage from the external power unit VP232-2. The fuse compartment of the VP232-2 power supply unit is located at the rear of the unit and contains two type T 4.0 A/H / 250 V fuses.

Pull out the power plug before changing fuses.

- Pull out the power plug.
- Pull out the fuse holder (Fig. 13/2) frontward.
- Remove the fuses from the fuse holder and insert new fuses.
- Reinsert the fuse holder into the fuse compartment (Fig. 13/1) until it clicks into place.
- Plug in the power plug.


Fig. 13 Changing the fuses in the Axio Observer 7 materials

#### 4.10.3 Changing the fuses in the ballast unit for the HBO 100

The fuse holder for the F1 and F2 fuses is located on the rear of the ballast unit. The fuse holder is integrated into the power socket and contains two T 2.0 A/H 250 V 5 x 20 mm fuses.

If the device fuses fail the cause must first of all be ascertained and technical problems properly remedied.

Pull out the power plug before changing fuses.

- Pull out the power plug.
- Remove the fuse holder (Fig. 14/1) from the fuse compartment (Fig. 14/2) by pulling it towards you.
- Replace the defective fuse.
- Reinsert the fuse holder into the fuse compartment until it clicks into place.
- Plug in the power plug.



Fig. 14 Changing fuses in the ballast unit

## 5 Shipping conditions

Pull out the power plug before transport.

- Switch off the Axio Observer (Fig. 15/1).
- Disconnect the device from the mains voltage by removing the power cable (Fig. 15/2) from the Axio Observer.
- Remove the transmitted light illuminator, microscope stage and binocular tube,
- Screw on the transport handles (Fig. 15/3) using a 4 mm Allen key.

Fig. 15 Transporting the Axio Observer

Full transportation descriptions are provided in the operating manual of the Axio Observer 431004-7144-001.

All Axio Observer models may only be transported individually across a few meters using the mounted transport handles (Fig. 15/3). For transporting over long distances the device must be partially dismantled and wrapped in the original packaging.

Following transportation, the device should be reassembled in the reverse order.

---

## 6 AMBIENT CONDITIONS

### Storage (in packaging)

Permissible ambient temperature.....	+5 °C to +40 °C
Permissible air humidity (without condensation).....	max. 75% at +35 °C

### Transport (in packaging):

Permissible ambient temperature.....	-40 °C to +70 °C
Permissible air humidity (without condensation).....	max. 75% at +35 °C

### Operation

Permissible ambient temperature.....	+10 °C to +35 °C, optimally 22 °C
Permissible relative air humidity .....	max. 65 % at 30 °C
Air pressure .....	800 hPa to 1060 hPa
Degree of pollution.....	2
Highest permitted altitude of use .....	max. 2000 m

## 7 TECHNICAL DATA

### Dimensions (width x depth x height)

Axio Observer 3 materials, 5 materials, 7 materials stand..... approx. 295 mm x 805 mm x max. 707 mm

### Weight

Axio Observer 3 materials..... approx. 27 kg

Axio Observer 5 materials ..... approx. 30 kg

Axio Observer 7 materials..... approx. 36 kg

### Operating data for Axio Observer 3 materials and 5 materials, manual with integrated power supply unit, and motorized Axio Observer 7 materials with external power supply unit VP232-2.

Operating area..... Enclosed rooms

Protection class ..... I

Ingress protection rating..... IP 20

Electrical safety .....DIN EN 61010-1 (IEC 61010-1)  
and CSA and UL regulations

Overvoltage category ..... II

Suppression of interference ..... acc. to EN 55011 Class B

Noise immunity..... acc. to DIN EN 61326-1

Line voltage (Axio Observer 3 materials and 5 materials) ... 100 V to 127 V and 200 V to 240 VAC  $\pm 10\%$

Line voltage of external power supply unit of Axio Observer 7 materials ..... 100 V to 240 VAC  $\pm 10\%$   
A change of the line voltage is not required!

Mains frequency ..... 50 Hz to 60 Hz

Power consumption of Axio Observer 3 materials and 5 materials, manual ..... max. 300 VA

Power consumption of external power supply unit of Axio Observer 7 materials, mot..... max. 190 VA

### Power supply unit (ballast unit) HBO 100

Operating area..... Enclosed rooms

Protection class ..... I

Ingress protection rating..... IP 20

Line voltage ..... 100 V to 240 VAC

Mains frequency ..... 50 Hz to 60 Hz

Power consumption if operating with HBO 100 ..... 155 VA

### Fuses in accordance with IEC 127

Microscope stand Axio Observer 3 materials and 5 materials, manual ..... T 5 A/H / 250 V, 5x20 mm

Power supply unit VP232-2 for Axio Observer ..... T 4.0 A/H / 250 V, 5x20 mm

HBO 100 power supply unit (ballast unit)..... T 2.0 A/H / 250 V, 5x20 mm

**Light sources**

HBO 50W/AC mercury vapor short-arc lamp

Output.....50 W

Average service life ..... 100 h

HBO 103 W/2 mercury vapor short-arc lamp..... 100 W

**Optical-mechanical data**Stand with stage focussing .....with coarse focusing drive approx. 2mm/rotation  
and fine focusing drive approx. 1/10 coarse/fine focus transmission ratio,  
Total travel approx. 10 mm, 13 mm also possible

Change of objective.....via 6-objective revolver

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