

Important Information: Windows 10 Upgrade Alert



All ZEISS microscope systems using the Windows 7 operating system are not compatible with the upcoming Windows 10 upgrade.

Upgrading your ZEISS system with Windows 10 could have negative effects, including complete system failure.

To avoid any disruptions with your ZEISS system, **we strongly urge you to not upgrade your system with Windows 10 when it becomes available on July 29, 2015.**

Please visit [this website](#) to:

- Determine if you have already signed-up for this Windows 10 upgrade
- Cancel your Windows 10 upgrade reservation
- Turn off automatic reminders from Microsoft to upgrade

If you have any questions or concerns, please contact our North American Product and Applications support team:

- support@zeiss.com
- 1-800-509-3905, M-F, 10AM-6PM Eastern Time

There may be a slightly longer response time. We will respond to your inquiry as quickly as possible.

Carl Zeiss Microscopy, LLC
One Zeiss Dr.
Thornwood, New York 10594
United States

Carl Zeiss Canada Ltd.
45 Valleybrook Dr.
Toronto, Ontario M3B 2S6
Canada

[Forward to colleague](#)